



Telelift - Laboratory sample
transport





Safe and hygienic laboratory sample transport with the Telelift UniCar

Optimizing logistics processes in laboratory buildings - safe, quick and reliable

In laboratory logistics, it is important that samples arrive reliably and quickly at their destination for analysis. At the same time, high hygiene standards apply and transports must be carried out gently so as not to contaminate the laboratory results.

How are laboratory samples transported by Telelift?

The Telelift UniCar System is a rail-bound transport system with three-dimensional routing and self-propelled transport trolleys that ensure automated, gentle and safe transport of laboratory samples. even across several floors or parts of buildings.

Safe and gentle transport with the Telelift UniCar.




Thanks to a swivel-mounted suspension of the container, the transported material is not moved when the trolley changes from vertical to horizontal or overhead transport, but always remains aligned in a balanced position. This means that even sensitive samples can be transported safely and gently.

The electronically monitored transport containers ensure hygienic and safe transport from the delivery point to the transport destination with complete transport documentation. TeleControl transport order management ensures fast processing of transport tasks by the shortest possible route, thus minimizing turnaround times in the laboratory.





Safe and hygienic laboratory sample transport with the Telelift UniCar



Learn more about the Telelift solution
for laboratory sample transport





Safe and hygienic laboratory sample transport with the Telelift UniCar



Maximum flexibility with just one system

The Telelift system offers maximum flexibility when it comes to how the specimens can be transported.

In our trolley with swiveling transport insert, laboratory samples can be transported both individually and as a bulk transport - packed, loose or in racks.

Loading with pre-commissioned laboratory samples is also possible without any problem. Thereby almost all common containers and laboratory sample carriers fit easily into the trolley.

The integrity of the laboratory samples is ensured during the entire transport. The gentle, automated transport and antimicrobial coating of the trolleys minimize the risk of cross-contamination and ensure a high standard of hygiene.

You too can benefit by using an automated transport solution from Telelift. The reduced waiting times and increased efficiency in sample transport ensure that your lab can always work at optimized performance.

The Telelift system takes over, e.g. the transport logistics for:



- Blood and urine samples
- Tissue samples (individual and bulk transport possible)
- Blood and plasma samples
- Sensitive samples that cannot be dispatched with other transport systems
- Consumables and supplies
- Laboratory documents



References laboratory transport

Telelift - A company with a success story

Telelift GmbH has been active in the hospital and laboratory sector since 1964. Our internal transport systems solve the challenges of transporting medical goods and samples in sometimes very confined spaces. The safety of the transported goods is always a top priority.

Telelift has installed more than 1,500 systems in more than 40 countries - including numerous hospital laboratories - and is certified according to DIN EN ISO 9001:2015.

Our goal is to work with you to find the optimal solution for your laboratory logistics or beyond.

An excerpt from our extensive reference list:

Benefits

- Safe, hygienic and gentle transportation of all laboratory goods and laboratory samples
- Antimicrobial coating of the trolleys
- Reduction of the risk of cross-contamination
- Optimization of handling times
- Seamless monitoring



Singapore General Hospital Laboratory

Distribution of internal and external laboratory samples from over 40 local outpatient clinics to several sub-labs.



Bristol Southmead Pathology

The Telelift system transports laboratory goods across several buildings and thus enables an efficient operation and a continuous workflow.



AUSL della Romagna

Transport of laboratory samples across several floors and with multiple stations.



Telelift GmbH

Frauenstraße 28
82216 Maisach / Germany

Phone +49 (0)8141 315 91-0

info@telelift-logistic.com
www.telelift-logistic.com



telelift
*Innovation for **Logistic** Solutions*

Mehr als 1.500 Installationen weltweit!
Mehr als 55 Jahre Erfahrung!

■ MADE
■ IN
■ GERMANY